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The screenshot shows the web interface for 'Individual SBII Calculation'. The header includes the email 'xinjie_hui@foxmail.com', the institution 'Chinese Academy of Medical Sciences and Peking Union Medical College, Beijing, China', and social media icons. The navigation menu has 'Home', 'Introduction', 'Tools', 'Download', and 'Contact us'. A search bar and a 'Batched SBII Calculation' button are also present. The main content area is titled 'Individual SBII Calculation' and is divided into two steps. Step 1, 'Upload EDF and XML Files', instructs users to wait after clicking 'UPLOAD FILES' and not to leave the page. It features two file selection buttons labeled '选择文件' (Select File) and an 'UPLOAD FILES' button. Step 2, 'Please select the Oxygen Saturation Channel and click 'SUBMIT' to generate SBII results', includes a dropdown menu for 'Select SaO2 Channel' and a 'SUBMIT' button.

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Step1: Upload EDF and XML Files

Please wait a few moments after clicking "UPLOAD FILES", duration may vary based on file size (< 120MB per upload).
Do not leave this page during the upload.

选择文件 test.edf 选择文件 test.xml **UPLOAD FILES**

Job ID: 8951703804741

Files have been uploaded successfully.

Available signal channels: ['SaO2', 'H.R.', 'EEG(sec)', 'ECG', 'EMG', 'EOG(L)', 'EOG(R)', 'EEG', 'AIRFLOW', 'THOR RES', 'ABDO RES', 'POSITION', 'LIGHT', 'NEW AIR', 'OX stat']

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Step2: Please select the Oxygen Saturation Channel and click 'SUBMIT' to generate SBII results

SaO2 ▼ **SUBMIT**

Job ID: 8951703804741

Database has been generated.

SBII result has been generated.


The final **SBII** and **pRED_3p** indices both need to be divided by Total Sleep Time (TST, hours), which can be readily obtained from the PSG report.

id	pRED_3p * TST	SBII * TST
test	0.742	72.316

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
The screenshot shows the 'Batched SBII Calculation' page. At the top, there is a navigation bar with links for Home, Introduction, Tools, Download, and Contact us. A search icon and a button for 'Individual SBII Calculation' are also present. The main heading is 'Batched SBII Calculation'. Below it, a message states: 'Please upload a Database (.db). The [User's Guide](#) below will provide a method for locally generating a database file.' There are two input fields: 'Organization or Hospital [Optional]' and 'Email Address [Optional]: SBII results will be sent to the email you provided as an attachment.' Below these is a file upload section with the text 'Upload a DB file:' and a button labeled '选择文件' (Select File), followed by the text '未选择任何文件' (No file selected). At the bottom, there are two buttons: 'RESET' and 'SUBMIT'.

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The screenshot shows the 'Thanks for using SBII' page. At the top right, there is a button labeled 'DOWNLOAD YOUR RESULTS'. Below this, the page displays job details: 'Job number: 2281703805101', 'Organization/Hospital:', and 'Email: xinjie_hui@foxmail.com'. A message indicates: 'Dec 29 07:11:52 webserver sendEmail[28188]: Email was sent successfully!'. Below this, a note states: 'The final SBII and pRED_3p indices both need to be divided by Total Sleep Time (TST, hours), which can be readily obtained from the PSG report.' At the bottom, there is a table with the following data:

id	pRED_3p * TST	SBII * TST
shhs1-200077	0.742	72.316
shhs1-200078	0.383	43.059
shhs1-200079	0.232	35.664

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Please upload a Database (.db). The [User's Guide](#) below will provide a method for locally generating a database file.

Organization or Hospital [Optional]

Email Address [Optional]: SBII results will be sent to the email you provided as an attachment.

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Thanks for using SBII. Your job 2281703805101 is finished! Please check the email attachment.

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out.txt
166.00B